

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-600: BS-600, BS-620;

2.BS-300: BS-300, BS-320;

9.BS-600M: BS-600M, BS-620M;

3.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-380: BS-380, BS-390;

11.BS-1000M: BS-1000M, BS-1100M;

5.BS-400: BS-400, BS-420;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-430: BS-430, BS-450,BS-460;

13.BS-2800M:BS-2600M, BS-2800M;

7.BS-480: BS-480, BS-490;

14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 061624003

: 2025-03-27

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

AFU

English	α -L-Fucosidase
Русский	α -L-фукозидаза
Português	α -L-Fucosidase
Español	α -L-fucosidasa
Italiano	α -L-Fucosidasi
Türkçe	α -L-Fukosidazdır

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
AFU	U/L	BS-200E ¹	/	/	/	/	/	BS-200E ¹	/	/	/	/	/
		BS-300 ²	/	/	/	/	/	BS-300 ²	/	/	/	/	/
		BS-330E ³	/	/	/	/	/	BS-330E ³	/	/	/	/	/
		BS-380 ⁴	/	/	/	/	/	BS-380 ⁴	/	/	/	/	/
		BS-400 ⁵	/	/	/	/	/	BS-400 ⁵	/	/	/	/	/
		BS-430 ⁶	/	/	/	/	/	BS-430 ⁶	/	/	/	/	/
		BS-480 ⁷	/	/	/	/	/	BS-480 ⁷	/	/	/	/	/
		BS-600 ⁸	/	/	/	/	/	BS-600 ⁸	/	/	/	/	/
		BS-600M ⁹	/	/	/	/	/	BS-600M ⁹	/	/	/	/	/
		BS-800 ¹⁰	/	/	/	/	/	BS-800 ¹⁰	/	/	/	/	/
		BS-1000M ¹¹	26.3	20.9	—	31.7	1.8	BS-1000M ¹¹	65.5	52.3	—	78.7	4.4
		BS-2000 ¹²	/	/	/	/	/	BS-2000 ¹²	/	/	/	/	/
		BS-2800M ¹³	/	/	/	/	/	BS-2800M ¹³	/	/	/	/	/

mindray
AFU Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT# 061624003

EXP 2025-03-27

