

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-600: BS-600, BS-620;

2.BS-300: BS-300, BS-320;

9.BS-600M: BS-600M, BS-620M;

3.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-380: BS-380, BS-390;

11.BS-1000M: BS-1000M, BS-1100M;

5.BS-400: BS-400, BS-420;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-430: BS-430, BS-450,BS-460;

13.BS-2800M:BS-2600M, BS-2800M;

7.BS-480: BS-480, BS-490;

14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 061625001

EXP: 2026-02-19

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

AFU

English	α -L-Fucosidase
Русский	α -L-фукозидаза
Português	α -L-Fucosidase
Español	α -L-fucosidasa
Italiano	α -L-Fucosidasi
Türkçe	α -L-Fukosidazdır

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
AFU	U/L	BS-200E ¹	23.6	18.8	—	28.4	1.6	BS-200E ¹	60.3	48.3	—	72.3	4.0
		BS-300 ²	27.2	21.8	—	32.6	1.8	BS-300 ²	68.1	54.3	—	81.9	4.6
		BS-330E ³	23.6	18.8	—	28.4	1.6	BS-330E ³	60.3	48.3	—	72.3	4.0
		BS-380 ⁴	24.8	19.7	—	29.9	1.7	BS-380 ⁴	63.7	50.8	—	76.6	4.3
		BS-400 ⁵	25.1	20.0	—	30.2	1.7	BS-400 ⁵	63.9	51.0	—	76.8	4.3
		BS-430 ⁶	26.2	20.8	—	31.6	1.8	BS-430 ⁶	66.6	53.1	—	80.1	4.5
		BS-480 ⁷	23.8	19.0	—	28.6	1.6	BS-480 ⁷	60.5	48.2	—	72.8	4.1
		BS-600 ⁸	24.8	19.7	—	29.9	1.7	BS-600 ⁸	62.8	50.2	—	75.4	4.2
		BS-600M ⁹	30.4	24.4	—	36.4	2.0	BS-600M ⁹	76.8	61.5	—	92.1	5.1
		BS-800 ¹⁰	23.8	19.0	—	28.6	1.6	BS-800 ¹⁰	60.5	48.2	—	72.8	4.1
		BS-1000M ¹¹	29.0	23.3	—	34.7	1.9	BS-1000M ¹¹	73.4	58.7	—	88.1	4.9
		BS-2000 ¹²	28.0	22.3	—	33.7	1.9	BS-2000 ¹²	71.0	56.6	—	85.4	4.8
		BS-2800M ¹³	29.0	23.3	—	34.7	1.9	BS-2800M ¹³	73.4	58.7	—	88.1	4.9

mindray
AFU Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 061625001  2026-02-19

