

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-200E: BS-200E, BS-220E,;

8.BS-600: BS-600, BS-620;

2.BS-300: BS-300, BS-320;

9.BS-600M: BS-600M, BS-620M;

3.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

10.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

4.BS-380: BS-380, BS-390;

11.BS-1000M: BS-1000M, BS-1100M;

5.BS-400: BS-400, BS-420;

12.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

6.BS-430: BS-430, BS-450,BS-460;

13.BS-2800M:BS-2600M, BS-2800M;

7.BS-480: BS-480, BS-490;

14.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT : 061625008

EXP: 2026-10-21

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)

AFU

English	α -L-Fucosidase
Русский	α -L-фукозидаза
Português	α -L-Fucosidase
Español	α -L-fucosidasa
Italiano	α -L-Fucosidasi
Türkçe	α -L-Fukosidazdır

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
AFU	U/L	BS-200E ¹	19.7	15.8	—	23.6	1.3	BS-200E ¹	54.0	43.2	—	64.8	3.6
		BS-300 ²	22.3	17.8	—	26.8	1.5	BS-300 ²	60.1	48.1	—	72.1	4.0
		BS-330E ³	19.7	15.8	—	23.6	1.3	BS-330E ³	54.0	43.2	—	64.8	3.6
		BS-380 ⁴	20.7	16.5	—	24.9	1.4	BS-380 ⁴	57.1	45.7	—	68.5	3.8
		BS-400 ⁵	20.6	16.4	—	24.8	1.4	BS-400 ⁵	56.4	45.0	—	67.8	3.8
		BS-430 ⁶	22.0	17.5	—	26.5	1.5	BS-430 ⁶	59.9	47.9	—	71.9	4.0
		BS-480 ⁷	19.9	16.0	—	23.8	1.3	BS-480 ⁷	54.2	43.4	—	65.0	3.6
		BS-600 ⁸	20.4	16.2	—	24.6	1.4	BS-600 ⁸	55.4	44.3	—	66.5	3.7
		BS-600M ⁹	24.8	19.7	—	29.9	1.7	BS-600M ⁹	67.7	54.2	—	81.2	4.5
		BS-800 ¹⁰	19.9	16.0	—	23.8	1.3	BS-800 ¹⁰	54.2	43.4	—	65.0	3.6
		BS-1000M ¹¹	23.8	19.0	—	28.6	1.6	BS-1000M ¹¹	64.8	51.9	—	77.7	4.3
		BS-2000 ¹²	23.3	18.5	—	28.1	1.6	BS-2000 ¹²	63.4	50.8	—	76.0	4.2
		BS-2800M ¹³	23.8	19.0	—	28.6	1.6	BS-2800M ¹³	64.8	51.9	—	77.7	4.3

mindray
AFU Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 061625008

2026-10-21

