

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.

Русский: Данные совпадают во всех группах.

Português: A dados de cada grupo é a mesma.

Español: la datos de cada grupo es la misma.

Italiano: la dati di ogni gruppo è la stessa.

Türkçe: her grubun veri aynıdır.

1.BS-120: BS-120, BS-130, BS-180, BS-190;

11.BS-430: BS-430, BS-450,BS-460;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

12.BS-480: BS-480, BS-490;

3.BS-200E: BS-200E, BS-220E,;

13.BS-600: BS-600, BS-620;

4.BS-240: BS-230, BS-240;

14.BS-600M: BS-600M, BS-620M;

5.BS-240E: BS240E, BS240Pro;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

6.BS-300: BS-300, BS-320;

16.BS-1000M: BS-1000M, BS-1100M;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

18.BS-2800M:BS-2600M, BS-2800M;

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");

9.BS-380: BS-380, BS-390;

19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

10.BS-400: BS-400, BS-420;

LOT : 061825003

EXP: 2026-08-16

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)
CO2				TBA	
English		Carbon Dioxide		Total Bile Acids	
Русский		Двуокись углерода		желчные кислоты	
Português		Dióxido de Carbono		Ácido Biliar Total	
Español		Dióxido de carbono		ácidos biliares totales	
Italiano		Anidride carbonica		acidi biliari totali	
Türkçe		Karbon Dioksit		Total Safra Asitleridir	

Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
CO2	mmol/L	BS-200 ²	12.3	9.9	—	14.7	0.8	BS-200 ²	17.4	13.8	—	21.0	1.2
		BS-200E ³	12.8	10.1	—	15.5	0.9	BS-200E ³	18.2	14.6	—	21.8	1.2
		BS-240 ⁴	12.1	9.7	—	14.5	0.8	BS-240 ⁴	17.4	13.8	—	21.0	1.2
		BS-240E ⁵	12.8	10.1	—	15.5	0.9	BS-240E ⁵	18.2	14.6	—	21.8	1.2
		BS-300 ⁶	11.8	9.4	—	14.2	0.8	BS-300 ⁶	17.5	13.9	—	21.1	1.2
		BS-330E ⁷	12.8	10.1	—	15.5	0.9	BS-330E ⁷	18.2	14.6	—	21.8	1.2
		BS-360E ⁸	12.8	10.1	—	15.5	0.9	BS-360E ⁸	18.2	14.6	—	21.8	1.2
		BS-380 ⁹	12.8	10.1	—	15.5	0.9	BS-380 ⁹	18.2	14.6	—	21.8	1.2
		BS-400 ¹⁰	12.8	10.1	—	15.5	0.9	BS-400 ¹⁰	18.2	14.6	—	21.8	1.2
		BS-430 ¹¹	12.8	10.1	—	15.5	0.9	BS-430 ¹¹	18.2	14.6	—	21.8	1.2
		BS-480 ¹²	12.8	10.1	—	15.5	0.9	BS-480 ¹²	18.2	14.6	—	21.8	1.2
		BS-600 ¹³	12.8	10.1	—	15.5	0.9	BS-600 ¹³	18.2	14.6	—	21.8	1.2
		BS-600M ¹⁴	12.8	10.1	—	15.5	0.9	BS-600M ¹⁴	19.1	15.2	—	23.0	1.3
		BS-800 ¹⁵	12.8	10.1	—	15.5	0.9	BS-800 ¹⁵	18.2	14.6	—	21.8	1.2

		BS-1000M¹⁶	12.8	10.1	—	15.5	0.9	BS-1000M¹⁶	18.6	15.0	—	22.2	1.2
		BS-2000¹⁷	12.2	9.8	—	14.6	0.8	BS-2000¹⁷	17.9	14.3	—	21.5	1.2
		BS-2800M¹⁸	12.8	10.1	—	15.5	0.9	BS-2800M¹⁸	18.6	15.0	—	22.2	1.2
Abbreviated name	Unit	Level: L/N					Level: H/P						
		Model	Assay Value	Range (Assay Value±3SD)		1 SD	Model	Assay Value	Range (Assay Value±3SD)		1 SD		
TBA	μmol/L	BS-120¹	20.6	16.4	—	24.8	1.4	BS-120¹	41.3	32.9	—	49.7	2.8
		BS-200²	20.6	16.4	—	24.8	1.4	BS-200²	41.3	32.9	—	49.7	2.8
		BS-200E³	20.7	16.5	—	24.9	1.4	BS-200E³	40.6	32.5	—	48.7	2.7
		BS-240⁴	20.6	16.4	—	24.8	1.4	BS-240⁴	41.3	32.9	—	49.7	2.8
		BS-240E⁵	21.3	17.1	—	25.5	1.4	BS-240E⁵	40.9	32.8	—	49.0	2.7
		BS-300⁶	21.4	17.2	—	25.6	1.4	BS-300⁶	41.6	33.2	—	50.0	2.8
		BS-330E⁷	20.7	16.5	—	24.9	1.4	BS-330E⁷	40.6	32.5	—	48.7	2.7
		BS-360E⁸	20.7	16.5	—	24.9	1.4	BS-360E⁸	40.6	32.5	—	48.7	2.7
		BS-380⁹	21.4	17.2	—	25.6	1.4	BS-380⁹	41.6	33.2	—	50.0	2.8
		BS-400¹⁰	21.4	17.2	—	25.6	1.4	BS-400¹⁰	41.6	33.2	—	50.0	2.8
		BS-430¹¹	20.9	16.7	—	25.1	1.4	BS-430¹¹	40.7	32.6	—	48.8	2.7
		BS-480¹²	21.4	17.2	—	25.6	1.4	BS-480¹²	41.6	33.2	—	50.0	2.8
		BS-600¹³	21.4	17.2	—	25.6	1.4	BS-600¹³	41.6	33.2	—	50.0	2.8
		BS-600M¹⁴	21.4	17.2	—	25.6	1.4	BS-600M¹⁴	41.6	33.2	—	50.0	2.8
		BS-800¹⁵	21.2	17.0	—	25.4	1.4	BS-800¹⁵	41.4	33.0	—	49.8	2.8
		BS-1000M¹⁶	21.4	17.2	—	25.6	1.4	BS-1000M¹⁶	41.6	33.2	—	50.0	2.8
		BS-2000¹⁷	21.2	17.0	—	25.4	1.4	BS-2000¹⁷	41.1	32.7	—	49.5	2.8
		BS-2800M¹⁸	21.4	17.2	—	25.6	1.4	BS-2800M¹⁸	41.6	33.2	—	50.0	2.8



CO2 and TBA Multi Control

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

LOT 061825003

 2026-08-16

