

Value sheet of Mindray BS Measurement System

Русский : Таблица результатов для системы BS компании Mindray
Português : Planilha de valores do Sistema de Medição BS da Mindray
Español : Hoja de valores del sistema de medición Mindray BS
Italiano : Scheda dei valori del sistema di misurazione BS di Mindray
Türkçe : Mindray BS Ölçüm Sistemi'nin değer sayfası



The data of each group is same.			Русский: Данные совпадают во всех группах.		
Português : A dados de cada grupo é a mesma.			Español : la datos de cada grupo es la misma.		
Italiano : la dati di ogni gruppo è la stessa.			Türkçe : her grubun veri aynıdır.		
1.BS-200: BS-200, BS-220, BS-330, BS-350;			10.BS-430: BS-430, BS-450,BS-460;		
2.BS-200E: BS-200E, BS-220E,;			11.BS-480: BS-480, BS-490;		
3.BS-240: BS-230, BS-240;			12.BS-600: BS-600, BS-620;		
4.BS-240E: BS240E, BS240Pro;			13.BS-600M: BS-600M, BS-620M;		
5.BS-300: BS-300, BS-320;			14.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;		
6.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")			15.BS-1000M: BS-1000M, BS-1100M;		
7.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with"W9-" and software version starts with "35.00");			16.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;		
8.BS-380: BS-380, BS-390;			17.BS-2800M:BS-2600M, BS-2800M;		
9.BS-400: BS-400, BS-420;			18.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.		

LOT : 145223007

 : 2025-08-17

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value±3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа ± 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise ±3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo ±3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione±3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri±3SD)

CHE	
English	Cholinesterase
Русский	Холинэстераза
Português	Colinesterase
Español	Colinesterasa
Italiano	Colinesterasi
Türkçe	Kolinesteraz

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)	1 SD	Model	Assay Value	Range(Assay Value±3SD)	1 SD
CHE	U/L	BS-200 ¹	/	/	/	BS-430 ¹⁰	/	/	/
		BS-200E ²	/	/	/	BS-480 ¹¹	/	/	/
		BS-240 ³	/	/	/	BS-600 ¹²	/	/	/
		BS-240E ⁴	/	/	/	BS-600M ¹³	/	/	/
		BS-300 ⁵	/	/	/	BS-800 ¹⁴	/	/	/
		BS-330E ⁶	/	/	/	BS-1000M ¹⁵	5915	4727 — 7103	396
		BS-360E ⁷	/	/	/	BS-2000 ¹⁶	/	/	/
		BS-380 ⁸	/	/	/	BS-2800M ¹⁷	/	/	/
		BS-400 ⁹	/	/	/				
		BS-200 ¹	/	/	/	BS-430 ¹⁰	/	/	/
	µkat/L	BS-200E ²	/	/	/	BS-480 ¹¹	/	/	/
		BS-240 ³	/	/	/	BS-600 ¹²	/	/	/
		BS-240E ⁴	/	/	/	BS-600M ¹³	/	/	/
		BS-300 ⁵	/	/	/	BS-800 ¹⁴	/	/	/

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value±3SD)				1 SD	Model	Assay Value	Range(Assay Value±3SD)				1 SD
		BS-330E ⁶	/	/	/	/	/	/	BS-1000M ¹⁵	98.8	78.9	—	118.6	6.6	
		BS-360E ⁷	/	/	/	/	/	/	BS-2000 ¹⁶	/	/	/	/	/	/
		BS-380 ⁸	/	/	/	/	/	/	BS-2800M ¹⁷	/	/	/	/	/	/
		BS-400 ⁹	/	/	/	/	/	/							