

Value sheet of Mindray BS Measurement System

Русский: Таблица результатов для системы BS компании Mindray

Português: Planilha de valores do Sistema de Medição BS da Mindray

Español: Hoja de valores del sistema de medición Mindray BS

Italiano: Scheda dei valori del sistema di misurazione BS di Mindray

Türkçe: Mindray BS Ölçüm Sistemi'nin değer sayfası

**The data of each group is same.****Português: A dados de cada grupo é a mesma.****Italiano: la dati di ogni gruppo è la stessa.**

1.BS-120: BS-120, BS-130, BS-180, BS-190;

2.BS-200: BS-200, BS-220, BS-330, BS-350;

3.BS-200E: BS-200E, BS-220E,;

4.BS-240: BS-230, BS-240;

5.BS-240E: BS240E, BS240Pro;

6.BS-300: BS-300, BS-320;

7.BS-330E:BS-330E(Serial Number starts with "XQ-"), BS-350E(Serial Number starts with "XS-")

8.BS-360E: BS-360E, BS-370E, BS-350S, BS-360S, BS-

330E(V35.00)(Serial Number starts with "W8-" and software version starts with "35.00"), BS-350E(V35.00) (Serial Number starts with "W9-" and software version starts with "35.00");

9.BS-380: BS-380, BS-390;

10.BS-400: BS-400, BS-420;

Русский: Данные совпадают во всех группах.**Español: la datos de cada grupo es la misma.****Türkçe: her grubun veri aynıdır.**

11.BS-430: BS-430, BS-450,BS-460;

12.BS-480: BS-480, BS-490;

13.BS-600: BS-600, BS-620;

14.BS-600M: BS-600M, BS-620M;

15.BS-800: BS-800, BS-820, BS-800M, BS-820M, BS-1800, BS-1800plus;

16.BS-1000M: BS-1000M, BS-1100M;

17.BS-2000: BS-2000, BS-2200, BS-2000M, BS-2200M;

18.BS-2800M:BS-2600M, BS-2800M;

19.For applicable models of the analyte, please refer to the reagent parameter sheet and instruction.

LOT: 145424003**: 2025-06-18**

English	Abbreviated name	Model	Unit	Assay Value	Range (Assay Value $\pm$ 3SD)
Русский	сокращенное наименование	модель	Прибор	Результат анализа	Диапазон(результат анализа $\pm$ 3CO)
Português	Nome abreviado	Modelo	Unidade	Valores da análise	Faixa(Valores da análise $\pm$ 3SD)
Español	nombre abreviado	modelo	Unidad	Valor de ensayo	Rango(Valor de ensayo $\pm$ 3SD)
Italiano	abbreviazione	modelli	Unità	Valori di dosaggio	Intervallo (valore di concentrazione $\pm$ 3 DS)
Türkçe	kısaltılmış ad	model	Ünite	Tayin Değeri	Aralık (Tayin Değeri $\pm$ 3SD)
FUN					
English			Fructosamine		
Русский			Фруктозамин		
Português			Frutosamina		
Español			Fructosamina		
Italiano			Fruttosamina		
Türkçe			Fruktozamin		

Abbreviated name	Unit	Model	Assay Value	Range(Assay Value $\pm$ 3SD)	1 SD	Model	Assay Value	Range(Assay Value $\pm$ 3SD)	1 SD
FUN	$\mu\text{mol/L}$	BS-120 ¹	459	390 — 528	23	BS-400 ¹⁰	457	388 — 526	23
		BS-200 ²	455	386 — 524	23	BS-430 ¹¹	471	399 — 543	24
		BS-200E ³	478	406 — 550	24	BS-480 ¹²	467	398 — 536	23
		BS-240 ⁴	455	386 — 524	23	BS-600 ¹³	464	395 — 533	23
		BS-240E ⁵	464	395 — 533	23	BS-600M ¹⁴	463	394 — 532	23
		BS-300 ⁶	467	398 — 536	23	BS-800 ¹⁵	467	398 — 536	23
		BS-330E ⁷	478	406 — 550	24	BS-1000M ¹⁶	456	387 — 525	23
		BS-360E ⁸	471	399 — 543	24	BS-2000 ¹⁷	472	400 — 544	24
		BS-380 ⁹	444	378 — 510	22	BS-2800M ¹⁸	456	387 — 525	23

mindray
FUN Control P

For use on: BS-1000M, BS-2000,
BS-2800M, BS-600M,
BS-800

 **LOT** 145424003

 **2025-06-18**

