

D-Dimer Control



Reference Value Table

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (low-value)	LOT 6224090116 2026-02-06	D-D	CX Series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	μg/mL	0.76	0.38 ~ 1.14
				ng/mL	760	380 ~ 1140

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (low-value)	LOT 6224090116 2026-02-06	D-D	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	μg/mL	0.85	0.43 ~ 1.28
				ng/mL	850	425 ~ 1275