

Reference Value Table

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (medium-value)	LOT 6324090116  2026-02-03	D-D	CX Series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	μg/mL	2.25	1.80 ~ 2.70
				ng/mL	2250	1800 ~ 2700

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (medium-value)	LOT 6324090116  2026-02-03	D-D	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	μg/mL	2.37	1.90 ~ 2.84
				ng/mL	2370	1896 ~ 2844