

Reference Value Table

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (medium-value)	LOT 6324110116  2026-02-26	D-D	CX Series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	μg/mL	2.40	1.92 ~ 2.88
				ng/mL	2400	1920 ~ 2880

Product Name	Batch code And Use-by date	Abbreviation/ Analyte	Applicable Instruments	Unit	Reference Value	Scope
D-Dimer control (medium-value)	LOT 6324110116  2026-02-26	D-D	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	μg/mL	2.45	1.96 ~ 2.94
				ng/mL	2450	1960 ~ 2940