

Fibrin(ogen) Degradation Products (FDP) Control



Reference Value Table

| Product Name            | Batch code And Use-by date                                                                                                                                      | Abbreviation/ Analyte | Applicable Instruments                                                                    | Unit  | Reference Value | Scope |       |       |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|-------------------------------------------------------------------------------------------|-------|-----------------|-------|-------|-------|
| FDP control (low-value) | <div><div>LOT</div><div>6624120116</div><div></div><div>2026-04-21</div></div> | FDP                   | CX Series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd | μg/mL | 7.66            | 3.83  | ~     | 11.49 |
|                         | mg/L                                                                                                                                                            |                       |                                                                                           | 7.66  | 3.83            | ~     | 11.49 |       |