

Reference Value Table

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation /Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	CX series Auto Coagulation Analyzers of Shenzhen Mindray Bio-medical Electronics Co., Ltd	PT	s	11.64	9.31 ~ 13.97
	8624090218		APTT	s	29.79	23.83 ~ 35.75
			TT	s	16.06	12.85 ~ 19.27
	2025-12-20		FIB	g/L	2.68	2.14 ~ 3.22
				mg/dL	268	214 ~ 322

Product Name	Batch code And Use-by date	Applicable Instruments	Abbreviation /Analyte	Unit	Reference Value	Scope
Coagulation control (level 1)	LOT	C3510 and ExC810 Auto Coagulation Analyzers of Beijing Mindray Medical Instrument Co., Ltd	PT	s	13.04	10.43 ~ 15.65
	8624090218		APTT	s	30.24	24.19 ~ 36.29
			TT	s	15.72	12.58 ~ 18.86
	2025-12-20		FIB	g/L	2.69	2.15 ~ 3.23
				mg/dL	269	215 ~ 323