

B30 CONTROL

FOR ARKRAY AUTO BLOOD COUNTER THINKA CB-1010



2025-09-02

ASSAY VALUES AND EXPECTED RANGES

Parameter	Low level		Normal level		High level	
	<div>LOT</div> BC20325L-1		<div>LOT</div> BC20325N-1		<div>LOT</div> BC20325H-1	
	Mean	Range	Mean	Range	Mean	Range
WBC $\times 10^9/L$	2.4	\pm 0.6	8.9	\pm 1.0	19.7	\pm 2.0
RBC $\times 10^{12}/L$	2.76	\pm 0.25	4.33	\pm 0.30	5.22	\pm 0.40
HGB g/L	83	\pm 5	133	\pm 7	167	\pm 9
HCT %	24.8	\pm 3.5	38.8	\pm 4.5	48.1	\pm 5.0
MCV fL	89.7	\pm 5.0	89.7	\pm 5.0	92.2	\pm 6.0
MCH pg	30.1	\pm 3.0	30.7	\pm 3.0	32.0	\pm 4.0
MCHC g/L	335	\pm 30	342	\pm 40	347	\pm 40
PLT $\times 10^9/L$	71	\pm 20	270	\pm 50	588	\pm 90

P/N:046-008372-00-1.0