

Hemoglobin A1c Control

mindray

Value sheet of Mindray Automated Glycohemoglobin Analyzer

LOT:202303071		🕒:2024-03-06					
Parameter	Method	Analyzer	Level	Unit		Assay Value	Range
HbA1c	HPLC	H50/H50P	CRL-1	IFCC	mmol/mol	39	33~45
				NGSP	%	5.7	5.2~6.2
			CRL-2	IFCC	mmol/mol	96	89~103
				NGSP	%	10.9	10.3~11.5