

Value sheet of Mindray Measurement System

LOT		060925001		2026-02-14				
Parameter	Method	Unit	Analyzer	Level	Values	Range		
C-reactive Protein (CRP)	Latex Immunoturbidimetric Method	mg/L	BC-5390 CRP、BC-5310 CRP BC-5100 CRP、BC-5180 CRP	Low	4.51	3.16	-	5.86
				High	51.67	36.17	-	67.17
C-reactive Protein (CRP)	Latex Immunoturbidimetric Method	mg/L	CRP-M100	Low	4.51	3.16	-	5.86
				High	51.67	36.17	-	67.17
C-reactive Protein (CRP)	Latex Particle-enhanced Immunonephelometric Method	mg/L	BC-7500[B] CRP、BC-7500[N] CRP BC-7500[R] CRP、BC-7500[NR] CRP BC-7500[B] CS、BC-7500[N] CS BC-7500[R] CS、BC-7500[NR] CS BC-7300[B] CRP、BC-7300[N] CRP BC-7300[R] CRP、BC-7300[NR] CRP	Low	5.47	3.83	-	7.11
				High	53.45	37.42	-	69.49