

Applicable Disinfecting Agents List

Disinfecting agent	Type	Manufacturer	ECG cable ECG lead set	Mobile monitoring ECG lead set	Temp Probe Temp adapter cable	IBP cable IBP adapter cable	C.O. cable assembly	NIBP cuff			PiCCO cable	SpO ₂ Sensor	
								CM12XX series/ CM17XX series	CM13XX series	CM1300C/ CM1300D/ CM1300E		512E/512M/512G/518B/518 BL/518BLH/513A/551B/51 2K/512RS	512F/512FLL/ 512FLH/512H
1-Propanol (50%)	Liquid	/	☒	☒	☒	☒	☒	☒	☒		☒	☒	☒
Alpet® D2 Surface Sanitizing Wipes	Wipes	BEST SANITIZERS Inc.	☒		☒	☒	☒	☒	☒		☒	☒	☒
CIDEX® OPA Solution	Liquid	Johnson & Johnson Medical GmbH	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Clorox Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Wipes	Clorox	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Clorox Healthcare® Bleach Germicidal Wipes	Wipes	Clorox	☒	☒	☒	☒	☒		☒	☒	☒	☒	☒
Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipe	Wipes	Clorox	☒	☒	☒	☒	☒		☒	☒	☒	☒	
Diversey Oxivir® TB Wipes	Wipes	Diversey Inc.	☒	☒	☒	☒	☒		☒		☒	☒	
Hydrogen peroxide (3%)	Liquid	/	☒	☒	☒	☒	☒		☒	☒	☒	☒	
Isopropanol (70%)	Liquid	/	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Metrex CaviCide 1™	Liquid	Metrex	☒	☒	☒	☒	☒		☒		☒	☒	
Metrex CaviWipes™	Wipes	Metrex	☒	☒	☒	☒	☒	☒	☒		☒	☒	☒
PDI Sani-Cloth® AF3 Germicidal Disposable Wipe	Wipes	PDI Inc.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
PDI Sani-Cloth® Bleach Germicidal Disposable Wipe	Wipes	PDI Inc.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
PDI Sani-Cloth® HB Germicidal Disposable Wipe	Wipes	PDI Inc.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
PDI Sani-Cloth® Plus Germicidal Disposable Cloth	Wipes	PDI Inc.	☒	☒	☒	☒	☒	☒	☒		☒	☒	☒
PDI Super Sani-Cloth® Germicidal Disposable Wipe	Wipes	PDI Inc.	☒	☒	☒	☒	☒	☒	☒		☒	☒	☒
Rely+On™ Virkon® Powder	Powder	DuPont	☒	☒	☒	☒	☒		☒	☒	☒	☒	
Sodium hypochlorite bleach (0.5%)	Liquid	/	☒	☒	☒	☒	☒		☒	☒	☒	☒	
VIRAGUARD Hospital Surface Disinfectants Towelettes	Wipes	VERIDIEN	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Virex® II 256	Liquid	Diversey Inc.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
Virex® TB	Liquid	Diversey Inc.	☒	☒	☒	☒	☒		☒	☒	☒	☒	
Mikrozid® sensitive liquid	Liquid	Schülke & Mayr GmbH	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
ANIOS Surf Premium	Liquid	ANIOS	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
ANIOS Surfa'safe Premium	Liquid	ANIOS	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
Clinell Universal Sanitizing Wipes	Wipes	GAMA Healthcare Ltd.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	
DIAN'ERKANG Surface Wipes	Wipes	Shanghai Likang Disinfectant Hi-Tech Co., Ltd.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
HEALTH ESSENCE Surface Disinfectant	Liquid	Beijing ChangJiangMai Medical Science Technology Co. Ltd.	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Glutaraldehyde (2%)	Liquid	/	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒
Ethanol (70%)	Liquid	/	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒	☒



消毒剂适用表

[illegible]