




## BeneHeart R12



Electrocardiograph



Innovations in R12 focus on improving healthcare quality, increasing operational efficiency and reducing user's costs with an affordable and compact solution

BeneHeart R12 delivers rapid and accurate ECG diagnosis to help you meet a variety of clinical demands



  Environment and Modernization

Paperless process delivers faster patient testing and reduces the use of costly ECG paper

- Full print preview image helps you decide if an immediate repeat of an ECG recording is needed, giving you a faster workflow than ever before
- Preview image for instantaneous on-screen review allows you to send data to the ECG management system without printing
- Review on screen helps you to retrieve previous ECG reports at any time



### Easy to use technology saves you time



- 8-inch high-resolution colour display
- Optional touchscreen



- User friendly workflow buttons
- Modern designed soft touch hotkeys
- Standard keyboard layout



Patient demographics can be retrieved from the worklist



Anatomically designed cable to minimise tangling and lead reversal



Optional barcode scanner for rapid input of accurate patient demographics

### Intelligent design for busy professionals



- Stable, clean, and accurate ECG waveforms ensure the quality of ECG data and aid faster processing
- Clear view of results helps you focus on diagnosis and care
- Save up to 800 ECG records internally, or use the USB flash drive for offline ECG storage and transmission
- Lithium ion battery powers greater than 3.5 hours of continuous operation or 400-page recording on one charge
- One of the lightest 12-lead electrocardiographs available – only 4.8 kg
- Optional hardy and compact trolley for easy mobility and storage

