


Maintain & update equipment remotely

The Mindray M-IoT Device Manager supports remote maintenance and management of your devices, maximising the service-life of crucial technology and protecting your investment. To help further maintain the reliability of your equipment, Mindray's technical team can provide comprehensive support.

- Export log files remotely for faster troubleshooting and Mindray support
- Remote software updates when idle help reduce maintenance interruptions
- Roll out configurations to selected devices simultaneously, providing maximum support for clinical procedures



M-IoT Device Manager

Data-driven device management





Maintenance through surveillance

Ensure the safety and effectiveness of all equipment at all times, with real-time device data on the M-IoT Device Manager.

Comprehensive equipment data

- Access device list and location records, for rapid equipment location and asset inventory
- Real-time failure monitoring and guidance for timely maintenance and reduced equipment downtime
- Battery status and replacement countdown, enhancing patient safety by limiting interruptions
- IP/MAC address for network access control to ensure cybersecurity



Category	Model	Serial Number	Status	Location	IP Address	MAC Address
Anesthesia	A7	SN_12345678	Working	OR Bed1	196.76.0.1	12:84:8F:8D:C2:78
Anesthesia	A7	SN_12345679	Working	OR Bed2	196.76.0.2	12:84:8F:8D:C2:79
Ventilator	SV800	SN_12345680	Working	ICU Bed1	196.76.0.3	12:84:8F:8D:C2:80
Ventilator	SV800	SN_12345681	Standby	ICU Bed2	196.76.0.4	12:84:8F:8D:C2:81
Monitor	N17	SN_12345682	Working	ICU Bed1	196.76.0.5	12:84:8F:8D:C2:82
Vital Sign	VS9	SN_12345683	Standby	GW ROOM1	196.76.0.6	12:84:8F:8D:C2:83

Designed to combat the everyday challenges faced by clinical engineers, M-IoT Device Manager supports the effective location and real-time availability of your equipment, helping reduce delays and maximise device up-time.

Achieving maximum usability

The M-IoT Device Manager obtains maintenance information from all connected equipment, providing statistical analysis on distribution and usage in various departments, helping protect your investment at all times.

- Streamline device allocation by monitoring usage analytics and flexibly assigning devices to different departments
- Confidently schedule and coordinate servicing thanks to greater device visibility
- Monitor the whole device life cycle to inform the purchase of replacement equipment

Complete device integration

M-IoT Device Manager supports Mindray and third-party equipment connected via BeneLink, providing a seamless device management solution for devices across your hospital.

Standard protocol extension

Simplify the management of all equipment with the integrated third-party asset management system through the SNMP interface.

