

# MX8 Cardiovascular Ultrasound System

## Transducer Family



### **CW2 (Non-imaging CW)**

Applications: Adult Cardiac,  
Pediatric Cardiac  
CW frequency: 2.0 MHz  
Biopsy guide: NA



### **P10-4s (Phased)**

Applications: Adult Cardiac,  
Pediatric Cardiac  
FOV: 90°  
Biopsy guide: NA



### **P7-3Ts (TEE)**

Applications: Adult Cardiac  
FOV: 90°  
Biopsy guide: NA



### **SP5-1Ns (Phased)**

Applications: Cardiac, Vascular  
FOV (max): 90°  
Biopsy guide: NGB-011, stainless  
steel, undetachable



### **L9-3s (Linear)**

Applications: Vascular  
FOV: 43 mm  
Convex radius: 20 mm  
Biopsy guide: NGB-034, multi-angle,  
reusable



### **L12-3RCs (Linear)**

Applications: Vascular  
FOV: 38 mm  
Biopsy guide: NGB-043, stainless  
steel, detachable; NGB-044, stainless  
steel, detachable, out-of-phase



### **P8-2s (Phased)**

Applications: Cardiac  
FOV: 90°  
Biopsy guide: NA



### **SC5-1Ns (Convex)**

Applications: Vascular  
FOV: 73 mm  
Convex radius: 60 mm  
Biopsy guide: NGB-022, stainless  
steel, detachable

#### **Mindray North America**

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2023 Mindray DS USA, Inc. Subject to change. 04/23 P/N: 0002-08-40625 Rev 1.0

# **mindray**

healthcare within reach