

TE5 Ultrasound System

Transducer Family



C5-2s (Convex)

Applications: Adult Abdominal, OB/GYN, Urology, Vascular, Nerve
FOV: 76°
Convex Radius: 51 mm
Needle guide: NGB-015, multi-angle, reusable



L20-5s (Linear)

Applications: MSK, Nerve, Small Parts, Superficial, Vascular
FOV: 28.5 mm
Needle guide: NA



C11-3s (Micro-Convex)

Applications: Abdominal, Cardiac, Vascular, Nerve
FOV: 100°
Convex radius: 15 mm
Needle guide: NGB-018, multi-angle, reusable



L12-3RCs (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 38 mm
Needle guide: NGB-043, multi-angle, reusable



L9-3s (Linear)

Applications: Vascular, MSK, Nerve, Small Parts
FOV: 43 mm
Needle guide: NGB-034, multi-angle, reusable



L12-4s (Linear)

Applications: Vascular, MSK, Nerve, Small Parts
FOV: 38 mm
Needle guide: NGB-007, multi-angle, reusable



L14-6Ns (Linear)

Applications: MSK, Nerve, Small Parts, Superficial, Vascular
FOV: 38°
Needle guide: NGB-007, multi-angle, reusable



L11-3VNs (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 38°
Needle guide: NA
eSpacialNavi™: Enabled



L7-3s (Linear)

Applications: Vascular, MSK, Nerve, Small Parts
FOV: 38 mm
Needle guide: NGB-007, multi-angle, reusable



7L4s (Linear)

Applications: MSK, Small Parts, Vascular
FOV: 38mm
Needle guide: NGB-007, multi-angle, reusable



L14-6s (Linear)

Applications: MSK, Nerve, Small Parts, Superficial, Vascular
FOV: 25mm
Needle guide: NGB-016, multi-angle, reusable



L14-5sp (Linear)

Applications: Intraoperative, MSK, Nerve, Superficial, Vascular
FOV: 25.3 mm
Needle guide: NA

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

©2024 Mindray DS USA, Inc. Subject to change. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. 01/24 P/N: 0002-08-40753 Rev 1.0

mindray
healthcare within reach